

Automating a Spatial Profiling Workflow To Explore The Effects Of Hypoxia In The Tumor Microenvironment In Head And Neck Cancer

Abstract

Background: Head and neck cancer (HNC) is a heterogeneous group of malignancies that arise from the mucosal surfaces of the upper aerodigestive tract. The tumor microenvironment (TME) of HNC is characterized by the presence of immune cells, stromal cells, and extracellular matrix components. A key feature of the TME is hypoxia, which promotes tumor growth, invasion, and metastasis by altering the expression of genes involved in angiogenesis, cell survival, and metabolism. Understanding the complex interplay between hypoxia and immune infiltrates in the TME of HNC is crucial for the development of novel therapeutic strategies for the treatment of this disease. Whole transcriptome analysis by digital spatial profiling is an excellent method of probing the TME, but assessing large cohorts can be time consuming. Automating a profiling workflow to reduce hands-on time and region of interest (ROI) selection bias will enable exploration of large cohorts to identify mechanisms of action, potential drug targets, and biomarkers.

Methods: We developed an optimized spatial multi-omic workflow to enable high-throughput spatial analysis on GeoMx® Digital Spatial Profiler (DSP) using the Whole Transcriptome Atlas (WTA) and immuno-fluorescent morphology markers: SYTO82 (nuclei), CAIX (hypoxia), pan-cytokeratin (epithelium), CD3 (T-cells). A.I.-based analysis (Oncotopix® Discovery) of serial section H&E images and GeoMx IF images was developed to identify ROIs for GeoMx collection. Immune hot and cold selection used leukocyte density; tumor/stromal interface selection used epithelial areas. Areas of illumination (AOI) were chosen using concentric CAIX expression gradients. Integrated analysis of digital images using Oncotopix Discovery and the whole transcriptome was done to assess the above TME compartments.

Results: Automated ROI placement based on tumor/stroma, hypoxia and immune infiltration and AI /Deep Learning based AOI segmentation reduced AOI selection time and improved accuracy of tissue compartment enrichment, especially between samples and tissue types. Automated development of hypoxia gradient-based AOI enabled a selection strategy not possible in the standard DSP software. Cell phenotyping using IF morphology scan was used to supervise cell deconvolution. DSP results correlated well with patient outcomes.

Conclusions: This work shows that ROI-based spatial analyses can be used to explore the effects of hypoxia levels on immune infiltration in HNC. Automated AI-based ROI selection provides a means of sampling relevant tumor subtypes based on hypoxia and immune infiltrate criteria in an unbiased, reproducible manner, and can provide a standardized, automated method for selecting ROIs and segmenting AOIs across a cohort of mixed tissue types and pathological subtype, improving throughput

Methods: GeoMx DSP together with Oncotopix Discovery

excluded

Inflamed

Immune desert

Inflamed (low)

Inflamed

Immune

25 low



Hypoxia^{low}/immune^{high} high

Hypoxia^{low}/immune^{high}

16 1 1 Hypoxia^{low}/immune^{high} high

11 0 1











Figure 1. 18 FFPE blocks from head and neck cohort were used for this study. Details of their HPV-positive and negative status, P16 expression, hypoxia and immune response status were obtained.

17 1 ¹ Hypoxia^{low}/immune^{high} high 2.5 80 high Inflamed

Mixed low 27.5 0 low

¹ Kelly Hunter , ¹ Joana Campos , ¹ Jack L. McMurray, ¹ Melanie Weigand, ¹ Chris Morse, ¹ Ketaki Hardikar, ² Jeni Caldara,² Regan Baird,² David Mason,³ Kyla Teplitz,² James Mansfield,⁴ Hisham Mehanna,⁴ Jill Brooks

¹ Propath, Hereford, UK, ² Visiopharm, Horsholm, Denmark, ³ Bruker, Seattle, WA, USA, ⁴ Institute of Head and Neck Studies and Education, University of Birmingham, Birmingham, UK

ROI & AOI selection







Figure 1: Using automated hotspot detection to select regions of interest (ROIs) for GeoMx. 1A shows a 4-color GeoMx immunofluorescence (IF) image. For this application, the four colors were Syto13 (nuclear marker, blue), cytokeratin (green), CA-9 (red) and CD3 (yellow). 1B shows the Visiopharm overlay of a Deep Learning APP trained for tumor segmentation (tumor = blue, stroma = green). A hotspot analysis was performed for CA-9 within tumor regions, revealing localized CA-9 expression around the edge of the tumor (1C). 1D shows 4 regions of interest that were chosen based on the CA-9 hotspot analysis. Each of these ROIs can be further analyzed for tumor/stroma (or other marker-based segmentations) to form the areas of illumination (AOIs) that will be collected on GeoMx.

Figure 2: Creating AOIs from CA-9 hotspots using two different methods. 2A shows a whole-slide 4-color IF images taken on a GeoMx of Sample 16. 2B shows a zoomed-in portion of the whole slide. For this sample, 7 circular hotspots were chosen to be ROIs. These can be seen superimposed on the color IF image. ROIs 17 and 3-7 were segmented using a paint-to-train deep-learning-based approach into tumor and stromal regions (see Figure 4). ROI 2 was created larger than the others, segmented into tumor and stroma, and then the tumor region was segmented into a very high CA-9 expression region and the rest. The high CA-9 region was expanded in 50 µm steps to form concentric regions around the CA-9 hotspot. 2F-2K shows an expanded view of this process on sample 14 ROI 17.

Figure 4: Creation of AOIs based on tumor/stroma tissue segmentation. This is the second method for the generation of AOIs. In this case, each circular ROI was seg--learning-based approach into tumor and stromal regions, each of which became an AOI for GeoMx data collection. 4A - 4D shows the sequence for ROI 2 of Sample 18. The circular ROI was segmented into tumor (blue) and stroma (green), and each AOI is shown in 4C (tumor) and 4D (stroma). Figures 4E -4H show the same sequence of steps for ROI 3 of Sample 18.. Figures 4I-L and 4M-P show the same sequence of steps for ROI 4 and ROI 3 of Sample 16.

Whole Transcriptomics Analysis

Al segmentation

Exclusion

Segmentation of cells expressing punctate staining patterns, such as RNAscope

Contouring around scattered cells or objects to collect local background signals

> Contouring around tissue regions by intensity or distance

Figure 3. Oncotopix AI segmentation unlocks new applications for GeoMx DSP. 3A-D shows a series of other uses of Oncotopix segmentation for GeoMx DSP. The flexibility of trainable deep-learning-based AI means that a wide range of AOIs can be created. **3A-3B** shows the segmentation of cells or tissues using deep learning without interference from tissue autofluorescence, which can be problematic for thresholding methods. **3C** shows finding cells based on punctate RNA stains; one can choose how many spots per cell are to be considered positive. **3D** shows using contouring to create a penumbra AOI around scattered cells or other objects, to enable the collection of local background signals. 3E shows using tissue segmentation or hotspot analysis to find a tissue type (in this case hypoxic tissue based on an IF marker) and then contours made around the central area.

Multiomics in Visiopharm

Multiomics—MASON

Summary

- Propath have demonstrated a deployable pipeline, integrating Oncotopix Discovery from Visiopharm and the GeoMx DSP software, leveraging Deep-Learning based image analysis upstream and downstream of the laboratory protocol.
- Manual adjustment of any ROIs or AOIs automatically generated in Oncotopix Discovery or GeoMx DSP prior to GeoMx DSP data collection is easy and powerful
- Annotation-trained AI/DL AOI segmentation is a significant improvement on traditional threshold-based methods, offering sensitivity, specificity and workflow improvements.

Fully automated ROI and AOI selection are possible for GeoMx DSP using On cotopix Discovery, which can greatly reduce the amount of technician and pathologist time required for large cohort studies

Use of AI to assist in ROI and AOI selection reduces selection bias and time and can increase cellular enrichment for profiling.

Hypoxia signature changes correlated with the contours generated based on hypoxia surrogate IHC marker, validating the spatial concentration gradient methodology

Increasing the cohort from 18 patients confirmed targets of interest from a 6 patient pi-

