

and Image Analysis (ESSB 2026)

Arutha Kulasinghe¹, Felicia Roland², Michelle Poulin², Ritu Mihani², Katherine Hales³, Daniel Winkowski³



Abstract

There remains a critical need for robust prognostic and predictive biomarkers to improve patient stratification for cancer immunotherapy across solid tumors. The tumour microenvironment (TME), particularly the tumour–stroma interface, is increasingly recognized as a potentially key determinant of therapeutic response, yet remains insufficiently characterized. In this study, we leveraged high-dimensional spatial proteomics combined with advanced image analytics to interrogate the cellular, functional, and metabolic architecture of tumour, stromal, and interface regions, with a focus on biologically active transition zones.

Using the Phenocycler-Fusion 2.0 platform with a custom immuno-metabolic panel targeting over 68 proteins, we profiled markers spanning tumour identity, stromal structure, immune cell subsets, functional states, and metabolic activity. This enabled high-resolution mapping of cell types, subtypes, and metabolic phenotypes across spatially defined compartments. Advanced image analysis via Visiopharm's pipeline facilitated accurate segmentation of tumour cores, stromal regions, and interface zones, while a Phenoplex™-guided workflow supported integrated cellular and metabolic characterization.

Our findings demonstrate that tumour and stromal core regions are both compositionally and metabolically distinct from their respective interface regions. While overall cellular proportions alone were not predictive of clinical outcomes, deeper analysis of functional and metabolic heterogeneity revealed important insights. Notably, we identified a tumour-specific metabolic signature associated with poorer clinical outcomes. In contrast, increased metabolic activity in immune cells localized to the tumour–stroma interface was associated with improved responses to immunotherapy, suggesting that spatially resolved immune fitness plays a central role in driving benefit.

Collectively, these results highlight the importance of integrating spatial proteomics with advanced image analytics to uncover biologically meaningful features of the TME. This approach provides a powerful framework for identifying novel biomarkers and improving our understanding of the mechanisms governing response to immunotherapy.

Results – Compositional differences between Primary and Metastatic TNBC




Figure 2. Representative H&E and mIF cores of TNBC TMA with annotated regions of interest (ROI).

[illegible]

Figure 3. Co-Localization Matrix shows cellular co-expression of biomarkers across the TME to validate accurate thresholding and increase confidence in phenotype outputs.




Figure 4. t-SNE maps of the cellular composition in the tumour (top) and stromal (bottom) regions between primary and metastatic TNBC

Results – TIL scoring




Figure 5. Immune cell scoring across each core for CD8, CD4, CD56, CD20, and across the tumor, stroma and margin regions.

Methods




Figure 1. H&E and representative mIF cores of TNBC tissue microarray (TMA), biomarker panel on PCF, and Phenoplex workflow.

Results – Metabolic scoring





Figure 6. Framework for metabolic marker scoring across the TME landscape.

Conclusions

This study highlights the framework for ultrahigh-plex discovery spatial proteomic profiling and image quantification in TNBC. Most notably, the modular analysis capability highlights how features of importance can be analysed (composition, TIL and metabolic scoring). Most notably, the metabolic scoring could be developed into a "H-score" enabling quantification of this signal to aid in functional cell phenotyping.