Deep Learning-Based, Fully Automated Analysis of Whole Slide Images Can Detect Invasive Breast Carcinoma and Count Ki-67 Easily and Precisely

/////VISIOPHARM®

Background:

- The basic treatment of hormone-receptor(HR) positive breast carcinoma(BC) is antihormonal therapy.
- Recently, cyclin-dependent kinases 4/6 inhibitor called abemaciclib, has been approved for high risk patients after surgery, together with **Ki-67** IHC MIB-1 pharmDx(Dako Omnis) for the patient selection.
- However, there remains concerns about reproducibility, and a development of automated image analysis is desired.

Conclusion:

- Here, we have created a fully automated software to analyze Ki-67, which includes **DL-based algorithm** to detect invasive BC nest and count a large amount of cells quickly.
- It has probed to have concordant result with experienced breast pathologists, although the prognostic value needs to be confirmed.
- As the software tend to have lower value than pathologists, we have to pay attention to the cut-off value.
- However, fully automated analysis can be performed by a technician or an engineer, and allow pathologists to focus on another task. This objective method would help clinicians and patients to select the optimal treatment

based on the tumor biology.


Design:

• Random **100 cases** were collected from a previous study [Histopathology. 2020;77(3):471-480]. (Figure 1)



Figure 1. Three serial sections were aligned by virtual triple staining method by Visiopharm Oncotopix software.

- Resected specimen of primary breast invasive carcinoma (HR positive and HER2 negative) had been stained with **Ki-67**(MIB-1; Dako, Autostainer Link 48; Agilent), and digitized by NanoZoomer-XR (Hamamatsu Photonics K.K.)
- Two board-certified pathologists evaluated Ki-67 following the recommendations from the International Working Group.(JNCI. 2011;103(22):1656-64. Figure 2)



- We also performed original scoring using a deep learning(DL)-based Ki-67(Breast, Al APP, research use only; Visiopharm) in a fully automated approach with no manual input or review.
- We compared the continuous output score with Spearman's correlation analysis, and calculated agreement with the clinically relevant cut-off of 20% for treatment with adjuvant abemaciclib.

Presented at USCAP 112th annual meeting, March 14, 2023 at New Orleans

Akira I. Hida ¹, Jeppe Thagaard ², Takahiro Watanabe ³, Thomas Ramsing ², Agnete Overgaard ², Yumi Oshiro ⁴ and Takuya Moriya ⁵ ¹ Matsuyama Shimin Hospital, ² Visiopharm A/S, ³ Hyogo College of Medicine, ⁴ Matsuyama Red Cross Hospital, ⁵ Kawasaki Medical School

Results:

- value was yielded by Path 1 and the software(AI-APP).
- **97.6%** respectively.





Figure 3. Each pair showed positive correlations, and AI-APP tend to yield lower value than pathologists(A-B). The threshold of 20% by Pathologist 1 was equivalent to 13.3% by AI-APP, which are indicated as orange dot line, and only four cases had a discordant result(B). When the same AI-APP was used by non-pathologists, the correlation was perfect(C).




 Median/mean value of each pathologist and automated analysis were 11.7/15.1(AIH, Path 1), 8.4/12.0(TW, Path 2) and 6.5/9.4 respectively. Every pair showed statistically positive correlations of Ki-67(Figure 3), and the highest

Time needed for pathologist to evaluate 100 cases ranged from 328(Path 2) to 562 minutes(Path 1). The software counted 385,985 cells on average within the full-face slides. When we set the threshold of Ki-67 at 20% by pathologist 1 and calculated 13% by software, the sensitivity and specificity of our new algorithm is 88.9% and

E-mail; ahida@matsuyama-shimin-hsp.or.jp / jth@visiopharm.com

Figure 4 & Table



	by	APP	
by Path 1	negative	positive	Total
Negative	3886	0	3886
weak +	316	1544	1860
Positive	14	5427	5441
Total	4216	6971	11187

Abstract #888

Figure 4. About 100 cells in each cases were annotated as positive/weak positive/ negative by the Pathologist 1(A: red circle; positive, yellow; weak positive, pink; negative), which were then compared with the result by the Al-APP(B: red dots; positive, blue dots; negative).

Out of total 7301 positive cells, the APP misses 330 cells (4.5% [2.9% among all 11187 cells]), and every single 3886 negative cells were judged as negative by the APP(Table).

A0054 Ki




Figure 5. Another pathologist assessed the ability of AI-APP about invasive cancer detection, as **good**(≥80% of preselected square of 1mm²)/fair(50-80%)/poor(<50%). Many cases were regarded as good(n=82: upper left), and 16/2 cases as fair/poor respectively(lower panel). Of note, APP seems to be superior to a pathologist on three "fair" cases, as AI-APP is more matched to the myoepithelial lining by p63(#A0054).