

Smriti Kala¹, Brenna O'Neill², Rasmus Lyngby², Rasmus Sorensen², Andrew Quong¹, Sam Lim¹, Clinton Hupple¹, Nina Lane¹, Jeppe Thagaard², Johan Dore Hansen², James Mansfield², Fabian Schneider², Michelle Macpherson¹ | 1) Standard BioTools Canada Inc., Markham, Ontario, Canada 2) Visiopharm A/S, Hørsholm, Denmark

Introduction

Imaging Mass Cytometry (IMC) is a powerful tool for the study of complex cellular interactions in the tumor microenvironment (TME) and in the discovery of biomarkers that can predict disease outcome or response to therapy. The Hyperion Imaging System (Standard BioTools) utilizes CyTOF[®] technology to simultaneously assess 40-plus protein markers at subcellular resolution without spectral overlap or background autofluorescence, thus providing unprecedented insight into the organization and function of the TME. However, despite the advances in staining and imaging methods, developments in analysis software have not kept pace, as we lack a complete, user-defined workflow in a single software package for analysis of high-plex imaging data. Here we present a comprehensive workflow using the Visiopharm[®] Phenoplex platform, designed specifically for high-plex IMC image analysis.

Methods and Materials

This study demonstrates a tissue phenotyping workflow in highly autofluorescent lung and colon cancer tissues using high-plex IMC, which offers the advantage of zero autofluorescence, and hence more reliable results. Serial sections of lung and colon cancer tissue microarrays (TMAs) were stained with a 26-marker panel comprised of structural, tumor, stroma, and immune cell markers (Table 1). The panel is an easily customized version of the Maxpar[®] Human Immuno-Oncology IMC Panel Kits. The IMC Cell Segmentation Kit (ICSK) was included in the panel for improved nucleus and plasma membrane demarcation. The data analysis pipeline uses Phenoplex software (Visiopharm) for easy, accurate, and quantifiable phenotyping. Analysis consisted of tissue segmentation, nuclear detection using a deep-learning algorithm pre-trained on IMC DNA channels, a threshold-based cellular phenotyping step, and spatial analyses.

Metal	Marker	Clone	Catalog Number
141Pr	αSMA	1A4	3141017D
143Nd	Vimentin	D21H3	3143027D
144Nd	CD14	EPR3653	3144025D
147Sm	MMP-13	OT12D8	Custom conjugation
148Nd	Pan-keratin	C11	3148020D
149Sm	CD11b	EPR1344	3149028D
151Eu	CD31	EPR3094	3151025D
152Sm	CD66b	BLR11H	91H033152
153Eu	CD44	IM7	3153029D
EPR20532-			
154Sm	EpCAM	222	91H024154
155Gd	FoxP3	PCH101	3155018D
156Gd	CD4	EPR6855	3156033D
158Gd	E-cadherin	24-E-10	3158029D
159Tb	CD68	KP1	3159035D
161Dy	CD20	HT	3161029D
162Dy	CD8a	C8/144B	3162035D
164Dy	iNOS	SP126	91H025164
166Er	CD45RA	HI100	3166031D
167Er	Granzyme B	EPR20129-217	3167021D
168Er	Ki-67	B56	3168022D
169Tm	Collagen 1	Polycolon	3169023D
170Er	CD3	Polycolon	3170019D
Cleaved caspase-3			
172Yb	5A1E	3172027D	
173Yb	CD45RO	UCHL1	3173016D
174Yb	HLA-DR	LIN3	3174025D
175Lu	Cyclin D1	SP4	Custom conjugation
195Pt	ICSK1		TIS-00001
196Pt	ICSK2		TIS-00001
198Pt	ICSK3		TIS-00001

Table 1. Custom Human Immuno-Oncology IMC Panel

Results

In this work, we have shown that analysis of IMC images from lung and colorectal cancer tissues can uncover tissue phenotypic signatures of the TME through the determination of immune cell types found in the vicinity of cancerous cells. Moreover, cell counts and tissue phenotypes were highly consistent across the serial sections, demonstrating the power of IMC in generating robust data.

Profiling the tumor microenvironment of colon and lung cancers




Figure 1. IMC staining of normal and cancerous colon and lung TMAs. Tumor tissue is identified by EpCAM+E-cadherin and pan-cytokeratin (pan-CK) staining. Mucinous adenocarcinoma of colon shows cells with co-localization of MMP-13 and granzyme B (yellow) surrounding the tumor cells, as do cells producing iNOS. Lung adenocarcinoma shows the presence of cyclin D1 and Ki-67 positive cells within the tumor. Scale bars are 200 μm. Enlarged versions of the selected areas are shown in the lower panels.

Cell segmentation using Phenoplex is consistent across serial sections




Figure 2. Tissue segmentation and cell segmentation using Phenoplex. To count cells in different tissue compartments, the tissue was differentiated into tumor (cyan) and stroma (blue). A 50 μm invasive margin (red) was depicted around the tumor region. Cell segmentation was based on pre-trained deep learning-based algorithms in the Phenoplex software, augmented with annotations from the IMC Cell Segmentation Kit. The enlarged area of the segmented image clearly shows the nucleus (blue) surrounded by the cellular margins (light grey).

Lung squamous cell carcinoma




Figure 3a. IMC staining of serial sections of lung squamous cell carcinoma. Scale bars are 200 μm.

Only 2.7% variability in cell counts across serial sections




Figure 3b. Serial sections of 10 different lung and colon cancer TMA cores were analyzed by the Phenoplex software. Cell segmentation showed equivalent number of cells between serial sections, with only 2.7% variability. This variability is not much greater than the 2% biological variability observed across serial sections in general. This result emphasizes the power of IMC in combination with analysis using Phenoplex software in generating highly consistent segmentation data using serial sections.

Lung large cell carcinoma-cellular clusters split by tissue compartments





Figure 4. Single-cell analysis of lung large cell carcinoma using Phenoplex can create t-SNE (t-distributed stochastic neighbor embedding) plots split by tissue compartments. The top left t-SNE (a) shows all lung cancer cells, and the top right t-SNE plots (b) show cellular clusters split according to the tissue compartments the cells were in. The bottom 12 t-SNE plots (c) show the cells with differential expression in the tumor, stromal and margin regions. For example, granzyme B and MMP-13 positive cells are in the tumor/margin but not stroma. Ki-67 positive cells are primarily in the tumor. T cell and other immune cell phenotypes vary by tissue region as well.

Conclusions

- This work demonstrates that even for highly autofluorescent tissues, IMC can generate high-quality data, consistent across serial sections.
- IMC data can be easily and accurately analyzed using Visiopharm's single software package Phenoplex, thus empowering IMC users to be confident in biological interpretation of high-dimensional proteomic data.
- This study showcases the capability of IMC technology combined with analysis using Phenoplex software for development of systematic digital profiling of the spatial TME.

Imaging Mass Cytometry workflow



Phenoplex: a new workflow for high-plex image analysis

